Type: Instructions for Use	Title: <b>SPECTRA Disposable Plastic and Metal Hub Epidural Needles - Instructions for Use</b>	Document No.: Rev. No.: Effective Date:	<b>SMD-IFU-0003 004 07 May 2021</b>
-------------------------------	---	---	---

SPECTRA Epidural Needle		
English	Instruction for Use	Epidural Needle

**Indications for Use:**

The **SPECTRA Epidural Needle** is intended for the epidural placement, directly or through an epidural catheter, of anesthetic agents to elicit regional anesthesia.

**Instructions for Use:**

- 1 Healthcare Practitioner must be trained / qualified in the use of the Epidural Needle for the intended procedure.
- 2 Read and follow all relevant hospital/medical institutional guidelines when placing an Epidural Needle for the subsequent medical procedure. Observe Universal Precautions on all patients.
- 3 Procedure should be carried out in an appropriate anesthetic environment.
- 4 Remove the needle Guard before use. Re-sheathing of needles is hazardous.
- 5 **NEVER PULL THE CATHETER BACK THROUGH THE EPIDURAL NEEDLE AS THIS MAY CAUSE THE CATHETER TO KINK OR SHEAR.**
- 6 Dispose of the SPECTRA Epidural Needle in an approved sharps container while observing universal precautions.
- 7 Report needle sticks immediately and follow relevant hospital/medical institutional protocol.
- 8 Report all adverse incidents involving these devices.

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0167-002-02, SMD-IFU-0003 Rev 004

Type: Instructions for Use	Title: SPECTRA Disposable Plastic and Metal Hub Epidural Needles - Instructions for Use	Document No.: Rev. No.: Effective Date:	SMD-IFU-0003 004 07 May 2021
-------------------------------	---	---	------------------------------------

SPECTRA- Epidural Needle	
	
Warnings:	Cautions:
<ul style="list-style-type: none"> <li>This product is for single use, <b>DO NOT REUSE - DO NOT RESTERILIZE</b>. Reuse and/or repackaging may create a risk of patient or user infection, compromise the structure integrity and/or essential material and design characteristics of the device, which may lead to device failure, and/or injury, illness or death of the patient.</li> <li>Do not use if package is opened or damaged.</li> <li>Do not alter this device or Straighten Bent Needle – Discard if Bent.</li> <li><b>DO NOT AUTOCLAVE</b></li> <li>This device may cause interference in proximity of a magnetic resonance imaging machine.</li> <li>Handle and discard in accordance with accepted medical practice and applicable regulations. After use this product is a potential biohazard.</li> <li><b>NEVER PULL CATHETER BACK THROUGH EPIDURAL NEEDLE</b> as this may cause the catheter to kink or shear.</li> <li>Re-sheathing of needles is hazardous.</li> </ul>	<ul style="list-style-type: none"> <li>Carefully read and follow all relevant hospital/medical institutional guidelines when placing an Epidural Needle for the subsequent medical procedure.</li> <li>Care and adherence to universal precautions must be taken to avoid accidental needle stick injuries.</li> <li>Only qualified healthcare practitioners with proper training should insert, manipulate and remove these devices.</li> <li>Follow standard infection control precautions as specified by the CDC (USA) or local equivalent.</li> <li>Federal (U.S.A.) Law restricts this device to sale by or on the order of a physician.</li> </ul>

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Contraindications – Absolute:	Contraindications – Relative:
<ul style="list-style-type: none"> <li>• Patient refusal</li> <li>• Coagulopathy</li> <li>• Local Skin infection at the injection site</li> <li>• Raised intracranial pressure</li> <li>• Local anaesthetic allergy</li> <li>• Uncooperative patient</li> <li>• Patients with significantly low Platelet Counts</li> <li>• Patients receiving Anticoagulation therapy</li> </ul>	<ul style="list-style-type: none"> <li>• Systemic infection</li> <li>• Fixed Cardiac output states</li> <li>• Uncorrected hypovolaemia</li> <li>• Spinal anatomical abnormalities</li> <li>• Pre-existing neurological disease</li> <li>• Prophylactic low-dose heparin therapy</li> <li>• Thrombocytopenia</li> </ul>

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Side Effects / Adverse Events:	Monitor in a Patient with an Epidural:
--------------------------------	--

- |   |  |
|---|--|
| <ul style="list-style-type: none"> <li>• Dural Puncture</li> <li>• Total Spinal Anesthesia</li> <li>• Nausea / vomiting</li> <li>• Topical anesthetic drug intoxication</li> <li>• Impairment of nerve or spine</li> <li>• Epidural hematoma / abscess</li> <li>• Infection</li> <li>• Headache</li> <li>• Cauda epidural symptoms</li> <li>• Transient neurological symptoms</li> <li>• Blood pressure fail</li> </ul> | <ul style="list-style-type: none"> <li>• Heart rate and blood pressure</li> <li>• Respiratory rate and oxygen saturation</li> <li>• Temperature</li> <li>• Sedation score</li> <li>• Pain score</li> <li>• Sensory level: height of the sensory block to cold</li> </ul> |
|---|--|

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