Type: Instructions for Use	Title: Spectra Guide Wire Introducer (GWI) Needle - Instructions for Use	Document No.: Rev. No.: Effective Date:	SMD-IFU-0002 004 28 Feb 2024
-------------------------------	---	---	---

1. PURPOSE

- 1.1. This document describes the requirements for Spectra Branded and Private Labeled Sterile Disposable Guide Wire Introducer (GWI) Needles Instructions for Use (IFU).

2. SCOPE

- 2.1. This document covers all packaged Spectra Branded and Private Labeled Sterile Disposable Guide Wire Introducer (GWI) Needles manufactured by Spectra Medical Devices, LLC (SMD).

3. APPROVED SUPPLIERS

- 3.1. Printing Solutions / MBM Packaging Labs, Inc.

4. ASSOCIATED DOCUMENTS

- 4.1. SMD-Spec-0033, Top Case Assembly Specification for Spectra Guide Wire Introducer (GWI) Needles
- 4.2. SMD-Spec-0022, Label Stock for Multi Language IFU's
- 4.3. SMD-Sec Spec-0004, Print Quality for SPECTRA Sterile Product Labels
- 4.4. SMD-Form-0097, Label Artwork Proofing and Approval Form
- 4.5. SMD-IP-0028, Inspection Procedure for Labels / Packaged Sterile Products

5. PROCEDURE

- 5.1. SPECTRA Label Stock part number per **SMD-Spec-0022**.
- 5.2. Printed graphics shall conform to SPECTRA specification **SMD-Sec Spec-0004**, Print Quality.
- 5.3. IFUs can be preprinted since it is common to all Spectra GWI Needles.
- 5.4. The IFU shall be printed in black ink (targeting PMS Cool Gray 11) 2-sides on 40 Lb. Opaque Uncoated Offset Paper Stock (12" x 13.3" sheet) to create a finished IFU Sheet.
- 5.5. The finished IFU shall be a 4.5" x 3" folded configuration.
- 5.6. The Printed Label must match the correct IFU configuration label shown in **EXHIBIT A**.

Type: Instructions for Use	Title: Spectra Guide Wire Introducer (GWI) Needle - Instructions for Use	Document No.: Rev. No.: Effective Date:	SMD-IFU-0002 004 28 Feb 2024
-------------------------------	---	---	------------------------------------

6. QUALITY REQUIREMENT

6.1. Refer to Label Artwork Proofing and Approval Form, **SMD-Form-0097**.

6.1.1. The “Verification of Label Text Requirements” section of this Form must be completed prior to release of a Purchase Order to print and fold the IFUs by the supplier. If accepted, the PO will be processed, and if rejected, the document will be corrected, and the process repeated.

6.1.2. After completing the Label Artwork Printable File, a copy will be inserted into the space provided in the “Label Artwork Review and Approval” section of the Form. The Form will then be routed for Review and Approval by a representative of at least two of the designated Management Functions and/or if required to the Customer. If accepted, the IFU proof file will be confirmed back to the supplier for release under their Doc Control System.

6.1.3. The original approved IFU proof file will be attached to the signed copy of form SMD-Form-0097 and filed.

6.1.4. In the event the Label Artwork is revised, this process must be repeated. In addition, the Supplier must be notified in writing to obsolete/remove the prior Approved Label Artwork and release the new Label Artwork approved by SMD.

6.2. **SMD-IP-0028**, Inspection Procedure for Packaged Sterile Products, will be utilized for Incoming Inspection of Sterile Packaged Product prior to Lot Release. A copy of the complete Record will be added to the DHR. A prior incoming inspection of the printed and folded IFUs shall be in place.

Type: Instructions for Use	Title: Spectra Guide Wire Introducer (GWI) Needle - Instructions for Use	Document No.: Rev. No.: Effective Date:	SMD-IFU-0002 004 28 Feb 2024
-------------------------------	---	---	------------------------------------

EXHIBIT A

INSTRUCTIONS FOR USE (IFU) LABEL ARTWORK

Type: Instructions for Use	Title: Spectra Guide Wire Introducer (GWI) Needle - Instructions for Use	Document No.: Rev. No.: Effective Date:	SMD-IFU-0002 004 28 Feb 2024
-------------------------------	---	---	------------------------------------

SPECTRA Guidewire Introducer Needle	
English	Instruction for Use Guidewire Introducer Needle
Indications for Use:	
The SPECTRA Guidewire Introducer Needle is indicated for percutaneous and placement of guidewires in vascular procedures.	
Instructions for Use:	
1	Remove the needle sheath.
2	Read and follow all instructions for use when placing a guidewire for the subsequent medical procedures.
3	Dispose of the SPECTRA Guidewire Introducer Needle in an approved sharps container while observing universal precautions.

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