



Guide Wire Introducer Needles




- Echogenic or Non-Echogenic Styles
- Available with or without Wing (base plate)
- Guaranteed Sharp Points (data available upon request)
- Special hub markings available for bevel indication
- Non-Coring "A" or "B" Bevels designed for anterior single wall arterial/seldinger percutaneous puncture
- Funneled, tapered luer hub and plug drawn tubing allows trouble-free guide wire insertion
- Available with clear hub or color coded by gauge
- Custom lengths available
- Calibrations available for the outer cannula



SIZE	GUIDE WIRE SIZE
18GAETW x 2 1/2" (7cm)	.038
18GAETW x 2 3/4" (9cm)	.038
19GAETW x 2 1/2" (7cm)	.032
19GAETW x 2 3/4" (9cm)	.032
19GAUTW x 2 1/2" (7cm)	.035
19GAUTW x 2 3/4" (9cm)	.035
20GAETW x 2 1/2" (7cm)	.025
20GAETW x 2 3/4" (9cm)	.025
21GATW x 2 1/2" (7cm)	.018
21GATW x 2 3/4" (9cm)	.018

Blunt Guidewire Introducer Needles

- Custom Lengths
- Multiple Gauges
- Round or Square Hub



Available with
Spectra Guide Glide Point

(978) 657-0889 • Fax: (978) 657-4339