



# *Potts Cournand & Modified Potts Cournand Style Needles*

- Sterile or Non-Sterile
- Funneled Tapered Luer Hub allows trouble-free Guide Wire Inserting (Modified Potts Cournand Only)
- Disposable Needle Set for Carotid Angiography as well as Direct Arterial Pressure Monitoring, Blood Sampling and Percutaneous Catherization



## **Potts Cournand Style**

18G (1.2mm) x 5.7cm  
.038 Guide Wire

with Obturator (3 part set)

18G (1.2mm) x 5.7cm  
.038 Guide Wire

## **Modified Potts Cournand Style**

18G (1.2mm) x 5.7cm  
.038 Guide Wire

with Obturator (3 part set)

18G (1.2mm) x 5.7cm  
.038 Guide Wire

**(978) 657-0889 • Fax: (978) 657-4339**

**260-F&H Fordham Road • Wilmington, MA 01887 • [www.SpectraMedical.com](http://www.SpectraMedical.com)**