



*Disposable*

# *Curved & Straight Radio Frequency Needles*

Designed for use with radiofrequency electrodes and generators in (P)RF procedures.

Cannula available in a variety of:

- Lengths
- Tip Sizes
- Gauges 18, 20, 21, 22
- Custom Gauges and Lengths available (see reverse side for more info)
- Safely compatible with all major brands of RF generators
- Machine Electro-Polished Tips for easy introduction

## **CURVED SHARP SPINAL TIP or BLUNT TIP**

- Curve allows for better cannula steerability with the green dot on the ergonomic hub indicating curve location & direction of the tip

## **STRAIGHT TIP**

- Ergonomic hub allows for easy advancement
- Quality insulation prevents RF leakage



**(978) 657-0889 • Fax: (978) 657-4339**

**260-F&H Fordham Road • Wilmington, MA 01887 • [www.SpectraMedical.com](http://www.SpectraMedical.com)**