



Spinal Needles

Quincke, Chiba & Whitacre Point Spinal Needles

**510(k)
CLEARED**

- 510(k) Cleared
- Available from 18 Ga to 27 Ga
- Lengths available from 4cm(1.5") to 30.5cm(12")
- Available with or without Introducers Needles
- Available with or without cm Markings
- Matched ground bevel
- Clear, Plastic Hubs with Color Coded Stylet Caps
- Bevel Orientation Indicators
- Luer Locking Hubs
- Protective Sheath
- Depth Guides available
- Depth Marks available



Long Length

Quincke & Chiba Point

Whitacre & All Other Pencil Point Styles

Aluminum Hub Spinal Introducer Needles

Call for a custom quote: 15 Ga to 22 Ga

- Female Luer
- Knurled Hub
- Guard for Protection
- Smooth Lead-In for needle passage / insertion
- Choice of Bevels
- Choice of Wall Thickness
- Several common lengths, e.g. 1-1/4" and 1-1/2"



(978) 657-0889 • Fax: (978) 657-4339

260-F&H Fordham Road • Wilmington, MA 01887 • www.SpectraMedical.com